


## YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 7mm; opening 5.8mm

110 gms; angled; fenestrated Ø 5mm



**SKU:** 244-541-07T

**Categories:** [Open Surgery](#), [Neurosurgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>MATERIAL</b>	<a href="#">Titanium</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.47.07
<b>SPECIFICATIONS</b>	<a href="#">angled</a> , <a href="#">temporary</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>MATERIAL</b>	<a href="#">Titanium</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.47.07

## SPECIFICATIONS

[angled](#), [temporary](#)