


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 10mm; opening 5.8mm

110 gms; angled; fenestrated Ø 5mm



SKU: 244-541-10T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.47.10
SPECIFICATIONS	angled , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.47.10

SPECIFICATIONS

[angled](#), [temporary](#)