


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 10mm; opening 5.8mm

110 gms; angled; fenestrated Ø 5mm




SKU: 244-543-10T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.47.20
SPECIFICATIONS	angled , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.47.20

SPECIFICATIONS

[angled](#), [temporary](#)