


## YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 11mm; opening 6.5mm

180 gms; curved



**SKU:** 244-403-11T

**Categories:** [Open Surgery](#), [Neurosurgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>MATERIAL</b>	<a href="#">Titanium</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.61.32
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">permanent</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>MATERIAL</b>	<a href="#">Titanium</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.61.32

## SPECIFICATIONS

[curved](#), [permanent](#)