


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 7mm; opening 7.8mm

110 gms; straight; fenestrated Ø 5mm



SKU: 244-545-07T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.47.27

SPECIFICATIONS [straight](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.47.27

SPECIFICATIONS

[straight](#), [temporary](#)