


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 10mm; opening 8.8mm

110 gms; straight; fenestrated Ø 5mm



SKU: 244-545-10T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.47.30

SPECIFICATIONS [straight](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.47.30

SPECIFICATIONS

[straight](#), [temporary](#)