


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 5mm; opening 5.8mm

110 gms; angled; fenestrated Ø 5mm



SKU: 244-547-05T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.47.35
SPECIFICATIONS	angled , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.47.35

SPECIFICATIONS

[angled](#), [temporary](#)