


## YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 7mm; opening 6.2mm

110 gms; angled; fenestrated Ø 5mm



**SKU:** 244-549-07T

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**MATERIAL** [Titanium](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.47.47

**SPECIFICATIONS** [angled](#), [temporary](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**MATERIAL** [Titanium](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.47.47

## SPECIFICATIONS

[angled](#), [temporary](#)