


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 7mm; opening 5.6mm

110 gms; angled; fenestrated Ø 5mm



SKU: 244-551-07T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.47.57
SPECIFICATIONS	angled , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.47.57

SPECIFICATIONS

[angled](#), [temporary](#)