


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 9mm; opening 7mm

180 gms; straight




SKU: 244-400-09T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT750T
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.60.10
SPECIFICATIONS	permanent , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT750T
MATERIAL	Titanium

AUTHORS

[YASARGIL](#)

Medicon

58.60.10

SPECIFICATIONS

[permanent](#), [straight](#)