


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 20mm; opening 13mm

180 gms: angled



SKU: 244-410-20T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.46
SPECIFICATIONS	angled , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.46

SPECIFICATIONS

[angled](#), [permanent](#)