


## YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 5mm; opening 5.6mm

200 gms; angled



**SKU:** 244-413-05T

**Categories:** [Open Surgery](#), [Neurosurgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>MATERIAL</b>	<a href="#">Titanium</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.61.55
<b>SPECIFICATIONS</b>	<a href="#">angled</a> , <a href="#">permanent</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>MATERIAL</b>	<a href="#">Titanium</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.61.55

## SPECIFICATIONS

[angled](#), [permanent](#)