


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 8mm; opening 5.4mm

200 gms; curved



SKU: 244-415-08T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.62
SPECIFICATIONS	curved , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.62

SPECIFICATIONS

[curved](#), [permanent](#)