


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 9mm; opening 7.2mm

200 gms; curved




SKU: 244-419-09T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT832T
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.68
SPECIFICATIONS	curved , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT832T
MATERIAL	Titanium

AUTHORS

[YASARGIL](#)

Medicon

58.61.68

SPECIFICATIONS

[curved](#), [permanent](#)