


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 15mm; opening 9.2mm

200 gms; straight



SKU: 244-400-15T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.60.16
SPECIFICATIONS	permanent , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.60.16

SPECIFICATIONS

[permanent](#), [straight](#)