


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 13mm; opening 7.5mm

200 gms; angled



SKU: 244-425-13T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.96
SPECIFICATIONS	angled , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.96

SPECIFICATIONS

[angled](#), [permanent](#)