


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 7mm; opening 7.4mm

150 gms; straight; fenestrated Ø 3.5mm




SKU: 244-430-07T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.59.50
SPECIFICATIONS	permanent , straight

Schließkraft
Closing Force



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.59.50

SPECIFICATIONS

[permanent](#), [straight](#)