


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 9mm; opening 8.4mm

150 gms; straight; fenestrated Ø 3.5mm



SKU: 244-430-09T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.59.51
SPECIFICATIONS	permanent , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.59.51

SPECIFICATIONS

[permanent](#), [straight](#)