


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 12mm; opening 9.6mm

180 gms; straight; fenestrated Ø 3.5mm



SKU: 244-430-12T
Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.59.52
SPECIFICATIONS	permanent , straight

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.59.52

SPECIFICATIONS

[permanent](#), [straight](#)