


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 7mm; opening 4.9mm

150 gms; angled; fenestrated \varnothing 3.5mm



SKU: 244-433-07T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.59.57
SPECIFICATIONS	angled , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.59.57

SPECIFICATIONS

[angled](#), [permanent](#)