


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 12mm; opening 10.3mm

180 gms; straight; fenestrated Ø 5mm



SKU: 244-436-12T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.59.62
SPECIFICATIONS	permanent , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.59.62

SPECIFICATIONS

[permanent](#), [straight](#)