


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 7mm; opening 7.2mm

150 gms; angled; fenestrated Ø 5mm



SKU: 244-437-07T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.59.64
SPECIFICATIONS	angled , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.59.64

SPECIFICATIONS

[angled](#), [permanent](#)