


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 10mm; opening 5.8mm

150 gms; angled; fenestrated Ø 5mm



SKU: 244-443-10T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.59.79
SPECIFICATIONS	angled , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.59.79

SPECIFICATIONS

[angled](#), [permanent](#)