


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 5mm; opening 7mm

150 gms; straight; fenestrated Ø 5mm



SKU: 244-445-05T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.59.81
SPECIFICATIONS	permanent , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.59.81

SPECIFICATIONS

[permanent](#), [straight](#)