


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 25mm; opening 13.3mm

200 gms; straight




SKU: 244-400-25T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.60.26
SPECIFICATIONS	permanent , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.60.26

SPECIFICATIONS

[permanent](#), [straight](#)