


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 7mm; opening 7.8mm

150 gms; straight; fenestrated Ø 5mm



SKU: 244-445-07T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.59.82
SPECIFICATIONS	permanent , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.59.82

SPECIFICATIONS

[permanent](#), [straight](#)