


## YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 7mm; opening 5.6mm

150 gms; angled; fenestrated Ø 5mm



**SKU:** 244-453-07T

**Categories:** [Open Surgery](#), [Neurosurgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)


**TYPE** [Vessel Clips](#)

**MATERIAL** [Titanium](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.59.94

**SPECIFICATIONS** [angled](#), [permanent](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**MATERIAL** [Titanium](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.59.94

## SPECIFICATIONS

[angled](#), [permanent](#)