


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 9mm; opening 7mm

90 gms; straight




SKU: 244-500-09T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT250T
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.60.60
SPECIFICATIONS	straight , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT250T
MATERIAL	Titanium

AUTHORS

[YASARGIL](#)

Medicon

58.60.60

SPECIFICATIONS

[straight](#), [temporary](#)