


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 9mm; opening 7mm

90 gms; straight



SKU: 244-500-09T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)


Aesculap FT250T

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.60.60

SPECIFICATIONS [straight](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FT250T

MATERIAL [Titanium](#)

AUTHORS

[YASARGIL](#)

Medicon

58.60.60

SPECIFICATIONS

[straight](#), [temporary](#)