


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 17mm; opening 10.6mm

90 gms; straight




SKU: 244-500-17T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.60.67
SPECIFICATIONS	straight , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.60.67

SPECIFICATIONS

[straight](#), [temporary](#)