

YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 20mm; opening 11.4mm

110 gms; straight



SKU: 244-500-20T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.60.71
SPECIFICATIONS	straight , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.60.71

SPECIFICATIONS

[straight](#), [temporary](#)