


## YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 9mm; opening 6.8mm

90 gms; curved



**SKU:** 244-501-09T

**Categories:** [Open Surgery](#), [Neurosurgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**Aesculap** FT252T

**MATERIAL** [Titanium](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.61.60

**SPECIFICATIONS** [curved](#), [temporary](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**Aesculap** FT252T

**MATERIAL** [Titanium](#)

## AUTHORS

[YASARGIL](#)

## Medicon

58.61.60

## SPECIFICATIONS

[curved](#), [temporary](#)