


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 9mm; opening 6.8mm

90 gms; curved



SKU: 244-501-09T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FT252T

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.61.60

SPECIFICATIONS [curved](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FT252T

MATERIAL [Titanium](#)

AUTHORS

[YASARGIL](#)

Medicon

58.61.60

SPECIFICATIONS

[curved](#), [temporary](#)