


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 9mm; opening 6.8mm

180 gms; curved



SKU: 244-401-09T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT752T
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.10
SPECIFICATIONS	curved , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT752T
MATERIAL	Titanium

AUTHORS

[YASARGIL](#)

Medicon

58.61.10

SPECIFICATIONS

[curved](#), [permanent](#)