


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 9mm; opening 6.1mm

110 gms; curved



SKU: 244-503-09T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.80
SPECIFICATIONS	curved , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.80

SPECIFICATIONS

[curved](#), [temporary](#)