


## YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 11mm; opening 6.2mm

110 gms; curved



**SKU:** 244-505-11T

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)


**TYPE** [Vessel Clips](#)

**MATERIAL** [Titanium](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.61.95

**SPECIFICATIONS** [curved](#), [temporary](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**MATERIAL** [Titanium](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.61.95

## SPECIFICATIONS

[curved](#), [temporary](#)