


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 7mm; opening 6.2mm

151 gms; straight



SKU: 244-400-07T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT740T
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.60.08
SPECIFICATIONS	permanent , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT740T
MATERIAL	Titanium

AUTHORS

[YASARGIL](#)

Medicon

58.60.08

SPECIFICATIONS

[permanent](#), [straight](#)