


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 7mm; opening 6.2mm

151 gms; straight



SKU: 244-400-07T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| Aesculap | FT740T |
| MATERIAL | Titanium |
| AUTHORS | YASARGIL |
| Medicon | 58.60.08 |
| SPECIFICATIONS | permanent , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|------------------------------|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| Aesculap | FT740T |
| MATERIAL | Titanium |

AUTHORS

[YASARGIL](#)

Medicon

58.60.08

SPECIFICATIONS

[permanent](#), [straight](#)