


YASARGIL VESSEL CLIP; STANDARD

permanent; 7mm; opening 6.2mm

150 gms; straight



SKU: 244-400-07

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE740

AUTHORS [YASARGIL](#)

SPECIFICATIONS [permanent](#), [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE740

AUTHORS [YASARGIL](#)

SPECIFICATIONS [permanent](#), [straight](#)